



OWEN MUMFORD
Making a World of Difference

HEAD OFFICE

INTERNATIONAL SALES & MARKETING

Owen Mumford Ltd
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : 01993 812021 F : 01993 813466
E : info@owenmumford.co.uk

UK / IRELAND SALES & MARKETING

Owen Mumford Ltd
Primsdown Industrial Estate
Worcester Road, Chipping Norton
Oxfordshire OX7 5XP, United Kingdom
T : 01608 645555 F : 01608 646546
E : info@owenmumford.co.uk

NORTH AMERICA

T : 800-421-6936
F : 770-977-2866
E : info@owenmumfordinc.com

FRANCE

T : 02 32 51 88 70
F : 02 32 21 96 95
E : contact@owenmumford-fr.com

GERMANY

T : 06026/97750
F : 06026/977529
E : info@owen-mumford.de

CHINA

T : 021 6032 3528
F : 021 6032 3782

www.owenmumford.com



Smartphone users can scan this code for more information about **AutoDrop** and **AutoSqueeze**.

OUR ACCREDITATIONS INCLUDE:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001



AT/1022-LT/MCC01128



AutoSqueeze[™]

AutoDrop[®]



OWEN MUMFORD



Simple eye drop delivery system

AutoDrop[®]

Eye Drop Guide




Using eye drops can be a challenge for many patients, young or old. **AutoDrop[®]** and **AutoSqueeze[™]** have been developed to make self administration as simple as possible – helping to improve patient compliance and reducing reliance on others to help with this task.

AutoDrop[®] and **AutoSqueeze[™]** are re-usable and provide a cost effective solution for simple eye drop application.

Unique design prevents blinking.

A pinhole draws the gaze away from the descending drop, whilst a small lip holds the lower eyelid down, this prevents blinking and ensures the eye drop enters the eye.



Compatible with most eye drop bottles.

AutoDrop[®] has been tested with hundreds of eye drops to check compatibility, for a full, up-to-date list please visit www.owenmumford.com or call us on **0800 032 0439**.



Delivers the drop in the right place every time.

The nozzle is positioned centrally over the eye ensuring that drops are correctly delivered on the first attempt.

Prevents contamination and cross infection.

The eye drop bottle is kept away from the eye at all times, preventing the tip from becoming contaminated.

AutoSqueeze[™]

Combine **AutoSqueeze[™]** with **AutoDrop[®]**, to make squeezing the bottle easier.



Many patients using eye drops also find it difficult to easily squeeze the eye drop bottle.

AutoSqueeze[™] clips around the bottle providing an ergonomic grip and extra leverage. It can be used in conjunction with **AutoDrop[®]** or on its own.



AutoDrop[®] and AutoSqueeze[™] have been designed to overcome common obstacles to patient compliance:

- Missing the eye with the drop.
- Blinking as the drop is delivered.
- Difficulty in squeezing the bottle to deliver the drop.

AutoDrop[®] is compatible with many eye drops on the market, including:

Antistin – Otrivine	Golden Eye – Golden Eye
Artelac – Artelac	Hypromellose – Tubilex
Azarga – Alcon	Liquifilm – Allergan
Betagan – Allergan	Lumigan – Allergan
Betnesol – UCB	Maxidrex- Alcon
Betoptic – Alcon	Maxitrol – Alcon
ChloromyctRedidrops – Goldshield	Murine – Prestige
Ciloxan – Alcon	Mydraciyl – Alcon
Combigan – Allergan	Nyogel – Novartis
Emadine – Alcon	Optrex – All Variants
Exocin – Allergan	Pred Forte – Allergan
Ganfort – Allergan	Timoptol – Teva
Genticin – Amdipharm	Xalatan – Pfizer

Tested as of 01/12/11

For the latest compatibility information, visit www.owenmumford.com or call us on **0800 032 0439** (Monday – Friday 8.30am – 5pm)